



The overall performance figures for standard weir-type HC4 Life Science line diaphragms includes:

- Recommended actuator operating pressure: 8.0 bar.G maximum
- Actuator temperature rating: -10° to 100°C
- Actuator autoclave temperature: 150°

EC Actuator (Standard 6 bar Spring) Performance Characteristics Rubber Diaphragms 100% ΔP Closure Condition



EC Actuator (Light 4 bar Spring) Performance Characteristics Rubber Diaphragms 100% ΔP Closure Condition



SC = Failsafe Spring Closing, SO = Failsafe Spring Opening, DA = Double Acting

www.cranechempharma.com

CRANE

CRANE ChemPharma Flow Solutions, Saunders® • 9860 Johnson Road
Montgomery, Texas USA 77316 • Tel.: +1 936 588 8360 • Fax: +1 936 588 8302

IMPORTANT: CRANE Co. and its subsidiaries are not responsible for the accuracy, compliance, and legality of material contained in this datasheet offered in print, on the company websites, or via any external links, or third party sites.